



Recombinant Human Calcium/calmodulin-dependent protein kinase type IV (CAMK4)

Product Code	CSB-MP614976HU
Abbreviation	CAMK4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16566
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLKVTVPSCS ASSCSSVTAS AAPGTASLVP DYWIDGSNRD ALSDFFEVES ELGRGATSIV YRCKQKGTQK PYALKVLKKT VDKKIVRTEI GVLLRLSHPN IIKLKEIFET PTEISLVLEL VTGGELFDRI VEKGYYSERD AADAVKQILE AVAYLHENG I VHRDLKPENL LYATPAPDAP LKIADFGLSK IVEHQVLMKT VCGTPGYCAP EILRGCA YGP EVD MWSV GII TYILLCGFEP FYDERGDQFM FRRILNCEYY FISPWWDEVS LNAKDLVRKL IVLDPKKRLT TFQALQHPWV TGKAANFVHM DTAQKKLQEF NARRKLKAAV KAVVASSRLG SASSSHGSIQ ESHKASRDPS PIQDGNEDMK AIPEGEKIQQ DGAQA AVKGA QAELMKVQAL EKVKGADINA EEAPKMV PKA VEDGIKVADL ELEEGLAE EK LKTVEEAAAP REGQGSSAVG FEVPQQDVIL PEY
Source	Mammalian cell
Target Names	CAMK4
Protein Names	Recommended name: Calcium/calmodulin-dependent protein kinase type IV Short name= CaMK IV EC= 2.7.11.17 Alternative name(s): CaM kinase-GR
Expression Region	1-473
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	The product of this gene belongs to the serine/threonine protein kinase family, and to the Ca(2+)/calmodulin-dependent protein kinase subfamily. This enzyme is a multifunctional serine/threonine protein kinase with limited tissue distribution, that has been implicated in transcriptional regulation in lymphocytes, neurons and male germ cells.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.