



Recombinant *Caenorhabditis elegans* mRNA-capping enzyme (cel-1), partial

Product Code	CSB-MP615037CXY
Abbreviation	cel-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17607
Product Type	Recombinant Protein
Immunogen Species	<i>Caenorhabditis elegans</i>
Purity	≥85% (SDS-PAGE)
Source	Mammalian cell
Target Names	cel-1
Protein Names	Recommended name: mRNA-capping enzyme Including the following 2 domains: Polynucleotide 5'-triphosphatase EC= 3.1.3.33Alternative name(s): mRNA 5'-triphosphatase Short name= TPase mRNA guanylyltransferase EC= 2.7.7.50Alt
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.