



# Recombinant Bovine Ubiquitin-conjugating enzyme E2 G2 (UBE2G2)

<b>Product Code</b>	CSB-YP615057BO
<b>Abbreviation</b>	UBE2G2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q17QG5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Bos taurus (Bovine)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGTALKRLM AEYKQLTLNP PEGIVAGPMN EENFFWEAL IMGPEdTCFE FGVFPAILS FPLDYPLSPPK MRFTCEMFHP NIYPDGRVCI SILHAPGDDP MGYESSAERW SPVQSVEKIL LSVVSMLAEP NDESGANVDA SKMWRDDREQ FYKVAKQIVQ KSLGL
<b>Source</b>	Yeast
<b>Target Names</b>	UBE2G2
<b>Protein Names</b>	Recommended name: Ubiquitin-conjugating enzyme E2 G2 EC= 6.3.2.19 Alternative name(s): Ubiquitin carrier protein G2 Ubiquitin-protein ligase G2
<b>Expression Region</b>	1-165
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein shares 100% sequence identity with the mouse counterpart. This gene is ubiquitously expressed, with high expression seen in adult muscle. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.