



# Recombinant Probable gamma-butyrobetaine dioxygenase (gbh-1)

<b>Product Code</b>	CSB-EP615182CXY-B
<b>Abbreviation</b>	gbh-1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q19000
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLSALLIRNI RNASKLASVA GPNSDRIVNV KWSDGKTGVF PLIWLRDTSP DPSTYTISPA MTARKLTMLE FDVEQNARKL WIDEDANCLK IEWESGVLSE FPSEWLKIRN PSDQEARRRR RKVYLFPEQT WGKAEIEGKL KKFSHEEFMK NEQVVHDFLQ AVCIDGIAVL KGAPQGVARGA VEAIGDRIGM IKRTHFGLVF EVSLKADASN MAYASNGGLP FHTDFPSLSH PPQLQMLHML QSAEEGGHSL FVDGFHVAEQ LRVEKPEIFK ILTTQSMEYI EEGYDVHEIN GKTIRFDYDM CARHKVIRLN DDGKVNKIQF GNAMRSWFYD CEP SKVQDVY RAMKTFTEYC YQPRNMLKFR LEDGDTVLWA NQRLLHTRDG FRNAPEKART LTGCYFDWDI VKSRVRFLRD KLSLEQNQPS A
<b>Source</b>	E.coli
<b>Target Names</b>	gbh-1
<b>Protein Names</b>	Recommended name: Probable gamma-butyrobetaine dioxygenase EC=1.14.11.1 Alternative name(s): Gamma-butyrobetaine hydroxylase Short name=Gamma-BBH Gamma-butyrobetaine,2-oxoglutarate dioxygenase
<b>Expression Region</b>	1-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.