



Recombinant Zinc finger protein unc-98 (unc-98)

Product Code	CSB-BP615185CXY
Abbreviation	unc-98
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19203
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MDDDFKEAR KERDDFEELM NACDLAKMSV KNNEMVHGLE TFGINGESSE NGNKEKPKEI MKVVAPTVEA YVGSSSAQTPTKSSGGALDG SDQQEVRQDG TSVQKDDNGF VFYKCRFCGL TFMFMNTLRA HERIHDVSRP YVCGKCGDSF EFACQLEYHA AQHSEIDGYK CECGRTFFSY TEMPLYHKHTD DPLELIGAPE TTTIKVSKKR VLPVSEQDLP QPAFVTEGYE PKHPLRVYND VRSKPYICEY CSKSYSDSRG LAYHMYSHRG EKYFNPRASR YMMGREGVGY TDSRSYYLFP RTSGYVTPRF
Source	Baculovirus
Target Names	unc-98
Protein Names	Recommended name: Zinc finger protein unc-98 Alternative name(s): Uncoordinated protein 98
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.