



Recombinant CRAL-TRIO domain-containing protein F28H7.8 (F28H7.8)

Product Code	CSB-BP615211CXY
Abbreviation	F28H7.8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q19895
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MAPNALSPVD IENMNDA AIE QVRLQVSDVI DPRYDTKW NM LRWLQSNDFN IPKTVHLLKK HLKWRKDRKL DEPESQSL LQ FSDARRKHAP IDIIGPQRKE DGDRLVVVDR AGRIDVSG LM KSVQPTEYLH EMFRSFEEIQ RRLMKMEAET GVQCYMHYIF DLEALNFDPT LLGVVNGPFR VSWQLVGQHY REFIDKFIVI NSPSYINVLW SALSPFIPEQ SKQRIVFAGS NWKEELLDIV DKECLPERYG GMIPDIQCLK PVDPIPKSLY WKLP AQYPTM DQLHKVSVSA SKHRMLIYKV DKPDTELLMY SHNENDITIT LYYSKNKNVS ENDL ELAVAP IPKCGLPAMD LFDYNCEYPG YYYIKLANEA SWLLPSTYRI IVIEKESGKE LEPLNLNEKW IKKGQKSKKK
Source	Baculovirus
Target Names	F28H7.8
Protein Names	Recommended name: CRAL-TRIO domain-containing protein F28H7.8
Expression Region	1-410
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.