



Recombinant Pipistrellus rusticus Forkhead box protein G1 (FOXG1)

Product Code	CSB-BP615223PBAF
Abbreviation	FOXG1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q1A1A4
Product Type	Recombinant Protein
Immunogen Species	Pipistrellus rusticus (Rusty pipistrelle)
Purity	>85% (SDS-PAGE)
Sequence	MLDMGDRKEV KMIPKSSFSI NSLVPEAVQN DNHHASHGHH NSHHPQHSHHH HHHHHHPPPP APQPPPPQ QPPPPQAPQP PQARGAPAAD DDKGPQQLLL PPPPPAAAL DGAKADGLGG KGEPGGGGAG ELAPVGPDEK EKGAGAGGEE KKGAGEGGKD GEGGKEGEKK NGKYEKPPFS YNALIMMAIR QSPEKRLTLN GIYEFIMKNF PYYRENKQGW QNSIRHNLSL NKCFVKVPRH YDDPGKGNYW MLDPSSDDVF IGGTTGKLRR RSTTSRAKLA FKRGARLTST GLTFMDRAGS LYWPMSPFLS LHHPRASSTL SYNGTTSAYP SHPMPYSSVL TQNSLGNHNS FSTANGLSVD RLVNGEIPYA THHLTAAALA ASVPCGLSVP CSGTYSLNPC SVNLLAGQTS YFFPHVPHPS MTSQSSTSMS ARAASSSTSP QAPSTLPCES LRPSLPSFTT GLSGGLSDYF THQNQGSSSN PLIH
Source	Baculovirus
Target Names	FOXG1
Protein Names	Recommended name: Forkhead box protein G1 Short name= FoxG1
Expression Region	1-484
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.