



Recombinant *Saccharomyces cerevisiae* Mediator of RNA polymerase II transcription subunit 2 (MED2)

Product Code	CSB-YP615463SVG
Abbreviation	MED2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12124
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVVQNSPVSS VHTANFSERG SNTRTMTYKN KLTVCFDDIL KVGAEEMMQQ QLKNVQLDSY LVNGFSQSQQ KLLKEKVKLF HGILDDLETS LSQSSSYLET LTALGKEKEK EREEAEKKRA EQENMRKVRE QEELKKRQEL EEASQQQQLQ QNSKEKNGLG LNFSTTAPAN TTDANGSKEN YQELGSLQSS SQTQLENANA ANNGAAFSP L TTRIQSQQA QPSDVMFNDL NSMDISMFSG LDSTGFDSTA FNATVDETKG FDDNDSGNNY NDINISSIEN NINNNINSTK NGKDNNNESN KNNNGDEKNK NNNEDNENNN NSSEKNNNNN NNNNNNDDN GNNNNNNSGN DNNNTTNDS NKNNSITTG NDNENIVNND LPTTVVSNPG DNPPPADNGE EYLTLNDFND LNIWSTTGD NGELDLSGFN I
Source	Yeast
Target Names	MED2
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 2 Alternative name(s): Mediator complex subunit 2
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.