



Recombinant *Saccharomyces cerevisiae* COP9 signalosome complex subunit 5 (RRI1)

Product Code	CSB-BP615508SVG
Abbreviation	RRI1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12468
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLSLNKTVKE LRQLLKERYT VEDELTESIA LSSMRFKPSQ EPEFHALSQS SLLKTKLKQQ SSTDIPSYTH VLISKLSCEK ITHYAVRGGN IEIMGILMGF TLKDNIVVMD CFNLPVVGTE TRVNAQLESY EYMVQYIDEM YNHNDGGDGR DYKGAKLNVV GWFHSHPGYD CWLSNIDIQT QDLNQRFDQP YVAIVVDPLK SLEDKILRMG AFRTIESKSD DNSATSYVEL ETIIFDSELN RALFETKLNL HCVIEDDESE QISLNRLIDS MKQYSYLMDS KNRVTRIKLA TTSERVSNEN KKNIDYQNRS TRSQFCLNTQ RGDSTETSSF GSMFSGDNTS DVDMEDRNLT EFDSTDTSLC INGEPsihvn RVERSSRSTD NFHNSKRMN SNQERCHDEG NDMLQRNVLE TDYARAKNRI LASKIKQYER LRFYKDTFTL
Source	Baculovirus
Target Names	RRI1
Protein Names	Recommended name: COP9 signalosome complex subunit 5 EC= 3.4.-.-
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.