



# Recombinant *Neurospora crassa* Chorismate synthase (aro-2)

<b>Product Code</b>	CSB-EP615522NHA-B
<b>Abbreviation</b>	aro-2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12640
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Neurospora crassa</i> (strain ATCC 24698 / 74-OR23-1A / CBS 708.71 / DSM 1257 / FGSC 987)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSTFGHYFRV TTYGESHCKS VGCIVDGVPP GMELTEDDIQ PQMTRRRPGQ SAITTPRDEK DRVIIQSGTE FGVTLGTPIG MLVMNEDQRP KDYGNKTM DI YPRPSHADWT YLEKYGVKAS SGGGRSSARE TIGRVAAGAI AEKYLKLAYG VEIVAFVSSV GSEHLFPPTA EHPSPSTNPE FLKLVNSITR ETVDSFLPVR CPDAEANKRM EDLITKFRDN HDSIGGTVTC VIRNVPSGLG EPAFDKLEAM LAHAML SIPA TKGFEVGS GF GGCEVPGSIH NDPFVSAENT EIPPSVAASG AARNGIPRPK LTTKTNFSGG IQGGISNGAP IYFRVGFKPA ATIGQEQT TA TYDGTSEGLV AAKGRHDPSV VPRAVPIVEA MAALVIMDAV LAQQARHTAK SLLPPLKQTI NSGKDTVGN G VSENVQESDL AQ
<b>Source</b>	E.coli
<b>Target Names</b>	aro-2
<b>Protein Names</b>	Recommended name: Chorismate synthase EC= 4.2.3.5 Alternative name(s): 5-enolpyruvylshikimate-3-phosphate phospholyase
<b>Expression Region</b>	1-432
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.