



Recombinant Human Ras-related protein Rab-40B (RAB40B)

Product Code	CSB-EP615534HU
Abbreviation	RAB40B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12829
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MSALGSPVRA YDFLLKFLLV GDSDVGKGEI LASLQDGAAE SPYGHPAGID YKTTTILLDG RRVKLQLWDT SGQGRFCTIF RSYSRGAQGV ILVYDIANRW SFDGIDRWIK EIDEHAPGVP KILVGNRLHL AFKRQVPTEQ AQAYAERLGV TFFEVSPLCN FNITESFTEL ARIVLLRHGM DRLWRPSKVL SLQDLCCRAV VSCTPVHLVD KLPLPIALRS HLKSFSMANG LNARMMHGGG YSLTTSSTHK RSSLRKVKLV RPPQSPPKNC TRNSCKIS
Source	E.coli
Target Names	RAB40B
Protein Names	Recommended name: Ras-related protein Rab-40B Alternative name(s): SOCS box-containing protein RAR Short name= Protein Rar
Expression Region	1-278
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	This protein has similarity to a yeast protein which suggests a role of the gene product in regulating secretory vesicles.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.