



# Recombinant Human Protein FAM53B (FAM53B)

<b>Product Code</b>	CSB-MP617903HU
<b>Abbreviation</b>	FAM53B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14153
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MVMVLSESLS TRGADSIACG TFSRELHTPK KMSQGPTLFS CGIMENDRWR DLDRKCPLQI DQPSTSIWEC LPEKDSSLWH REAVTACAVT SLIKDLSISD HNGNPSAPPS KRQCRSLFS DEMSSCRTSW RPLGSKVWTP VEKRRCYSGG SVQRYSNFGS TMQRSSSFSL PSRANVLSSP CDQAGLHHRF GGQPCQGVPG SAPCGQAGDT WSPDLHPVGG GRLDLQRSLS CSHEQFSFVE YCPPSANSTP ASTPELARRS SGLSRRSRQP CVLNDKKVGV KRRRPEEVQE QRPSLDLAKM AQNCQTFSSL SCLSAGTEDC GPQSPFARHV SNTRAWTALL SASGPGGRTP AGTPVPEPLP PSFDDHLACQ EDLSCEESDS CALDEDCGRR AEPAAAWRDR GAPGNSLCSL DGELDIEQIE KN
<b>Source</b>	Mammalian cell
<b>Target Names</b>	FAM53B
<b>Protein Names</b>	Recommended name: Protein FAM53B
<b>Expression Region</b>	1-422
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.