



Recombinant Human Protein FAM53B (FAM53B)

Product Code	CSB-MP617903HU
Abbreviation	FAM53B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14153
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MVMVLSESLS TRGADSIACG TFSRELHTPK KMSQGPTLFS CGIMENDRWR DLDRKCP LQI DQPSTSIWEC LPEKDSSLWH REAVTACAVT SLIKDLSISD HNGNPSAPPS KRQCRSLSFS DEMSSCRTSW RPLGSKVWTP VEKRRCYSGG SVQRYSNNGFS TMQRSSSFSL PSRANVLSSP CDQAGLHHRF GGQPCQGVPG SAPCGQAGDT WSPDLHPVGG GRLDLQRSLS CSHEQFSFVE YCPPSANSTP ASTPELARRS SGLSRRSRQP CVLNDKKVGV KRRRPEEVQE QRPSLDLAKM AQNCQTFSSL SCLSAGTEDC GPQSPFARHV SNTRAWTALL SASGPGGRTP AGTPVPEPLP PSFDDHLACQ EDLSCEESDS CALDEDGRR AEPAAAWRDR GAPGNSLCSL DGELDIEQIE KN
Source	Mammalian cell
Target Names	FAM53B
Protein Names	Recommended name: Protein FAM53B
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.