



Recombinant Human Four and a half LIM domains protein 2 (FHL2)

Product Code	CSB-MP617906HU
Abbreviation	FHL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14192
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MTERFDCHHC NESLFGKKYI LREESPYCVV CFETLFANTC EECGKPIGCD CKDLSYKDRH WHEACFHCSQ CRNSLVDKPF AAKEDQLLCT DCYSNEYSSK CQECKKTIMP GTRKMEYKGS SWHETCFICH RCQQPIGTKS FIPKDNQNF VPCYEKQHAM QCVQCKKPIT TGGVTYREQP WHKECFVCTA CRKQLSGQRF TARDDFAYCL NCFCDLYAKK CAGCTNPISG LGGTKYISFE ERQWHNDCFN CKKCSLSLVG RGFLTERDDI LCPDCGKDI
Source	Mammalian cell
Target Names	FHL2
Protein Names	Recommended name: Four and a half LIM domains protein 2 Short name= FHL-2 Alternative name(s): LIM domain protein DRAL Skeletal muscle LIM- protein 3 Short name= SLIM-3
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.