



Recombinant Human Hyaluronan-binding protein 2 (HABP2), partial

Product Code	CSB-MP617919HU
Abbreviation	HABP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14520
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	FSLMSLL ESLDPDWTPD QYDYSYEDYN QEENTSSTLT HAENPDWYYT EDQADPCQPN PCEHGGDCLV HGSTFTCSCL APFSGNKCQK VQNTCKDNPC GRGQCLITQS PPYYRCVCKH PYTGPCSCSQV VPVCRPNPCQ NGATCSRHKR RSKFTCACPD QFKGKFCEIG SDDCYVGDGY SYRGKMNRTV NQHACLYWNS HLLLQENYNM FMEDAETHGI GEHNFCRNPD ADEKPWCFIK VTNDKVKWEY CDVSACSAQD VAYPEESPTE PSTKLPGFDS CGKTEIAERK IKR
Source	Mammalian cell
Target Names	HABP2
Protein Names	Recommended name: Hyaluronan-binding protein 2 EC= 3.4.21.-Alternative name(s): Factor VII-activating protease Factor seven-activating protease Short name= FSAP Hepatocyte growth factor activator-like protein Plasma hyaluronan-b
Expression Region	24-313
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.