



Recombinant Human N-alpha-acetyltransferase 30 (NAA30)

Product Code	CSB-MP617933HU
Abbreviation	NAA30
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q147X3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEVPPGPSS LLPPPAPPAP AAVEPRCPFP AGAALACCSE DEEDDEEHEG GGSRSPAGGE SATVAAKGHP CLRCPPQEQE QQQLNGLISP ELRHLRAAAS LKSKVLSVAE VAATTATPDG GPRATATKGA GVHSGERPPH SLSSNARTAV PSPVEAAAAS DPAAARNGLA EGTEQEEEEEE DEQVRLSSS LTADCSLRSP SGREVEPGED RTIRYVRYES ELQMPDIMRL ITKDLSEPYI IYTRYFIHN WPQLCFLAMV GEECVGAIVC KLDMHKKMFR RGYIAMLAVD SKYRRNGIGT NLVKKAIYAM VEGDCDEVVL ETEITNKSAL KLYENLGFVR DKRLFRYYLN GVDALRLKLW LR
Source	Mammalian cell
Target Names	NAA30
Protein Names	Recommended name: N-alpha-acetyltransferase 30 EC= 2.3.1.88 Alternative name(s): N-acetyltransferase 12 N-acetyltransferase MAK3 homolog NatC catalytic subunit
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.