



# Recombinant Human ETS-related transcription factor Elf-2 (ELF2)

<b>Product Code</b>	CSB-EP618015HU-B
<b>Abbreviation</b>	ELF2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15723
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MTSAVVDSGG TILELSSNGV ENQEESKVS EYPAVIVEPV PSARLEQGYA AQLVYDDET YMMQDVAEEQ EVETENVETV EASVHSSNAH CTDKTIEAAE ALLHMESPTC LRDSRSPVEV FVPPCVSTPE FIHAAMRPDV ITETVVEVST EESEPMDTSP IPTSPDSHEP MKKKKVGKRP KTQQSPISNG SPELGIKKKP REGKGNTTYL WEFLDLLQD KNTCPRIKW TQREKGIFKL VDSKAVSKLW GKHKNKPD MN YETMGRALRY YYQRGILAKV EGQRLVYQFK DMPKNIVVID DDKSETCNED LAGTTDEKSL ERVSLSAESL LKAASSVRSG KNSSPINC SR AEKGVARVVN ITSPGHDASS RSPTTTASVS ATAAPRTVRV AMQVPVVM TS LGQKISTVAV QSVNAGAPLI TSTSPTTATS PKVVIQTIPT VMPASTENG D KITMQPAKII TIPATQLAQC QLQTKSNLTG SGSINIVGTP LAVRALTPVS IAHGTPVMRL SMPTQQASGQ TPRRVISAVI KGPEVKSEAV AKKQEHDVKT LQLVEEKPAD GNKTVTHVVV VSAPSAIALP VTMKTEGLVT CEK
<b>Source</b>	E.coli
<b>Target Names</b>	ELF2
<b>Protein Names</b>	Recommended name: ETS-related transcription factor Elf-2 Alternative name(s): E74-like factor 2 New ETS-related factor
<b>Expression Region</b>	1-593
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.