



Recombinant Human ETS-related transcription factor Elf-2 (ELF2)

Product Code	CSB-YP618015HU
Abbreviation	ELF2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15723
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTSAVVDSGG TILELSSNGV ENQESEKVS EYPAVIVEPV PSARLEQGYA AQLVYDDET YMMQDVAAEQ EVETENVETV EASVHSSNAH CTDKTIEAAE ALLHMESPTC LRDSRSPVEV FVPPCVSTPE FIHAAMRPDV ITETVVEVST EESEPMDTSP IPTSPDSHEP MKKKKVG RKP KTQQSPISNG SPELGIIKKP REGKGNTTYL WEFLDLLQD KNTCPRIKW TQREKGIFKL VDSKAVSKLW GKHKNKPD MN YETMGRALRY YYQRGILAKV EGQRLVYQFK DMPKNIVVID DDKSETCNE D LAGTTDEKSL ERVSLSAESL LKAASSVRSG KNSSPINC SR AEKGVARVVN ITSPGHDASS RSPTTTASVS ATAAPRTVRV AMQVPVVM TS LGQKISTVAV QSVNAGAPLI TSTSPTTATS PKVVIQTIPT VMPASTENG D KITMQPAKII TIPATQLAQC QLQTKSNLTG SGSINIVGTP LAVRALTPVS IAHGTPVMRL SMPTQQASGQ TPRRVISAVI KGPEVKSEAV AKKQEHDVKT LQLVEEKPAD GNKTVTHVVV VSAPSAIALP VTMKTEGLVT CEK
Source	Yeast
Target Names	ELF2
Protein Names	Recommended name: ETS-related transcription factor Elf-2 Alternative name(s): E74-like factor 2 New ETS-related factor
Expression Region	1-593
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.