



# Recombinant Human Cytochrome P450 1B1 (CYP1B1)

<b>Product Code</b>	CSB-EP618091HU-B
<b>Abbreviation</b>	CYP1B1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16678
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	<p>MGTSLSPNDP WPLNPLSIQQ TLLLLLSVL ATVHVGQRLL RQRRRQLRSA  PPGPFAPWLI GNAAAVGQAA HLSFARLARR YGDVFQIRLG SCPIVVLNGE  RAIHQALVQQ GSAFADRFAP ASFRVVSQGR SMAFGHYSEH  WKVQRRRAHS MMRNFFTRQP RSRQVLEGHV LSEARELVAL  LVRGSADGAF LDPRPLTVVA VANVMSAVCF GCRYSHDDPE FRELLSHNEE  FGRTVGAGSL VDVMPWLQYF PNPVRTVFRE FEQLNRNFSN FILDKFLRHC  ESLRPGAAPR DMMDAFILSA EKKAAGDSHG GGARLDLENV PATITDIFGA  SQDTLSTALQ WLLLLFTRYD DVQTRVQAEI DQVVGRDRLP CMGDQPNLPY  VLAFLYEAMR FSSFVPVTIP HATTANTSVL GYHIPKDTVV FVNQWSVNHD  PLKWPNPENF DPARFLDKDG LINKDLTSRV MIFSVGKRRRC IGEELSKMQL  FLFISILAHQ CDFRANPNEP AKMNFSYGLT IKPKSFKVNV TLRESMELLD  SAVQNLQAKE TCQ</p>
<b>Source</b>	E.coli
<b>Target Names</b>	CYP1B1
<b>Protein Names</b>	Recommended name: Cytochrome P450 1B1 EC= 1.14.14.1 Alternative name(s): CYP1B1
<b>Expression Region</b>	1-543
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The enzyme encoded by this gene localizes to the endoplasmic reticulum and metabolizes procarcinogens such as polycyclic aromatic hydrocarbons and 17beta-estradiol. Mutations in this gene have been



associated with primary congenital glaucoma; therefore it is thought that the enzyme also metabolizes a signaling molecule involved in eye development, possibly a steroid.

---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.