



Recombinant SUMO-activating enzyme subunit aos-1 (aos-1)

Product Code	CSB-MP618154CXY
Abbreviation	aos-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17820
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MEVSKAEQAI YDRQIRLWGM EAQNKIRNSK VLIIGGKQLG AEVAKTLLSLA GVDEMHLVDH RLVDTTEEIGM NFLYDASVDN SKMTKWAASY NFLYNLNRNV KLFIVEEDVL SKNDSEIEEY LTKFTLVVVL DESYERTAKV NNICRKHHR FISGAIYGIWY GYAFFDFDGH AYLKAKSPD CLNEEESETG KTSTVVTVDE EFVLETFSYP SFVETLNSDF TAKKIVRKCK RIVPTSIFLV KSMLRASSEN KLTGVTENDI EKLIPIWNEE VAAGNHTIDM QPVQPDRFDH LFGPNFGPTA ACVGGVIGQE AIKSISEGKN PLRNLFYITG FESTGFMCNF PPV
Source	Mammalian cell
Target Names	aos-1
Protein Names	Recommended name: SUMO-activating enzyme subunit aos-1
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.