



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YLL056C (YLL056C)

Product Code	CSB-MP618575SVG
Abbreviation	YLL056C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12177
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKVFITGASG FIGSAVLSEL ISSGHEVVGL ARSDEAAAKI KSIDPAAKIL RGDLKDLEIL KKGATESDGV IHLGFVHDFK NFEQCCEIDR QATVAMLESL KGSNKPFLYT NGTSLSRPNK VANEQDGIDE DSKILRAVTE QVALSYKDKG VSARIVRLPF SVHGKGDKAF VPILMNIKA AGKSGYVGQG TNAWAAVHRL DTAPLFRLVL EKGKTGQVYH CVGEQGIPFK DIARVIGEIL NVPVASIPVD DAESHFGFLT CFVTRDGPVS SEGTRKELGW QPQQIGLLED IRANYSLN
Source	Mammalian cell
Target Names	YLL056C
Protein Names	Recommended name: Uncharacterized protein YLL056C
Expression Region	1-298
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.