



Recombinant *Saccharomyces cerevisiae* Transcription initiation factor TFIID subunit 3 (TAF3)

Product Code	CSB-EP618589SVG
Abbreviation	TAF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12297
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MTTNNDFYFA LLRISILQLL KAQGFDRARP SLVDVMTDLY AKFLSLLASE VSSIAQARCD QDDTIALQDI TLALENLIV KPTNVLDVYD ENSELSSSRG MEKFKDWCIY STQLTDARIT ALPTVELLQS EEKESDPLSA IPDYLNQLLQ NKGAKQKLET KNRKTELIED LINNGLDDW IKLVIAQRRI NMIERASKKE SQNVPALPHI AGYKSSILSR HHHTTITNED RMPSAMTPRD EDALTEIQEN PFVTSKLPIM RKENRLENIT LSFEDDEELES LGEVEGPNQK SQENNNEESF KENNKSLTES PHGDDRDISM FQFDSNVDTK WAEQEDMDST FQRRTSLDYG GYF
Source	E.coli
Target Names	TAF3
Protein Names	Recommended name: Transcription initiation factor TFIID subunit 3 Alternative name(s): TAFII-47 Short name= TAFII47 TBP-associated factor 3 TBP-associated factor 47 kDa
Expression Region	1-353
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.