



Recombinant *Saccharomyces cerevisiae* Chromatin modification-related protein EAF3 (EAF3)

Product Code	CSB-MP618600SVG
Abbreviation	EAF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12432
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVDLEQEFAL GGRCLAFHGP LMYEAKILKI WDPSSKMYTS IPNDKPGGSS QATKEIKPQK LGEDESIPEE IINGKCFPIH YQGWKSSWDE WVGYDRIRAY NEENIAMKKR LANEAKEAKK SLLEQQKKKK LSTSLGGPSN GGKRKGDSRS NASISKSTSQ SFLTSSVSGR KSGRSSANSL HPGSSLRSSS DQNGNDDRRR SSSLSPNMLH HIAGYPTPKI SLQIPIKLS VLVDWEYVT KDKKICRLPA DVTVEMVLNK YEHEVSQELE SPGSQSQLSE YCAGLKLYFD KCLGNMLLYR LERLQYDELL KKSSKDQKPL VPIRIYGAIH LLRLISVLPE LISSTTMDLQ SCQLLIKQTE DFLVWLLMHV DEYFNDKDPN RSDDALYVNT SSQYEGVALG M
Source	Mammalian cell
Target Names	EAF3
Protein Names	Recommended name: Chromatin modification-related protein EAF3 Alternative name(s): ESA1-associated factor 3
Expression Region	1-401
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.