



Recombinant *Saccharomyces cerevisiae* Vacuolar-sorting protein SNF8 (SNF8)

Product Code	CSB-BP618605SVG
Abbreviation	SNF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12483
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MKQFGLAAFD ELKDGKYNDV NKTILEKQSV ELRDQLMVFQ ERLVEFAKKH NSELQASPEF RSKFMHMCSS IGIDPLSLFD RDKHLFTVND FYEVCCLKVI EICRQTKDMN GGVISFQELE KVHFRKLNVG LDDLEKSIDM LKSLECFEIF QIRGKKFLRS VPNELTSDQT KILEICSI LG YSSISLLKAN LGWEAVRSKS ALDEMVANGL LWIDYQGGAE ALYWDPSPWIT RQL
Source	Baculovirus
Target Names	SNF8
Protein Names	Recommended name: Vacuolar-sorting protein SNF8 Alternative name(s): ESCRT-II complex subunit VPS22 Vacuolar protein-sorting-associated protein 22
Expression Region	1-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.