



Recombinant *Saccharomyces cerevisiae* Actin-like protein ARP6 (ARP6)

Product Code	CSB-BP618609SVG
Abbreviation	ARP6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12509
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	METPPIVIDN GSYEIKFGPS TNKKPFRALN ALAKDKFGTS YLSNHIKNIK DISSITFRRP HELGQLTLWE LESCIWDYCL FNPSEFDGFD LKEGKGHHLV ASESCMTLPE LSKHADQVIF EEEYFDSLK SPVAVFVPFT KSYKGEMRTI SGKDEDIDIV RGNSDSTNST SSESKNAQDS GSDYHDFQLV IDSGFNCTWI IPVLKGIPYY KAVKKLDIGG RFLTGLLKET LSF RHYNMMD ETILVNNIKE QCLFVSPVSY FDSFKTKDKH ALEYVLPDFQ TSFLGYVRNP RKENVPLPED AQIITLDEL FTIPETFFHP EISQITKPGI VEAILESLSM LPEIVRPLMV GNIVCTGGNF NLPNFAQRLA AELQRQLPTD WTCHVSVPEG DCALFGWEVM SQFAKTDSYR KARVTREEYY EHGPDWCTKH RFGYQNW I
Source	Baculovirus
Target Names	ARP6
Protein Names	Recommended name: Actin-like protein ARP6
Expression Region	1-438
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.