



Recombinant Human Nuclear inhibitor of protein phosphatase 1 (PPP1R8)

Product Code	CSB-MP618641HU
Abbreviation	PPP1R8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12972
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGGEDDELKGLLGLPEEETELDNLTEFNHTAHNKRISTLTIEEGNLDIQRPKRKR KNSRVTFSEDDEIINPEDVDPSVGRFRNMVQTAVVPVKKKRVEGPGSLGLEES GSRRMQNFAFSGGLYGGLPPTHSEAGSQPHGIHGTALIGGLPMPYPNLAPDV DLTPVVPSAVNMNPAPNPAVYNPEAVNEPKKKKYAKEAWPGKKPTPSLLILEN GTHMASSDACHECKSMIAAAG
Source	Mammalian cell
Target Names	PPP1R8
Protein Names	Recommended name: Nuclear inhibitor of protein phosphatase 1 Short name= NIPP-1 Alternative name(s): Protein phosphatase 1 regulatory inhibitor subunit 8 Including the following 1 domains: Activator of RNA decay EC= 3.1.4.-A
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length of Isoform 2
Target Details	This gene, through alternative splicing, encodes three different isoforms. Two of the protein isoforms encoded by this gene are specific inhibitors of type 1 serine/threonine protein phosphatases and can bind but not cleave RNA. The third protein isoform lacks the phosphatase inhibitory function but is a single-strand endoribonuclease comparable to RNase E of E. coli. This isoform requires magnesium for its function and cleaves specific sites in A+U-rich regions of RNA.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.