



Recombinant Human Cleavage stimulation factor subunit 3 (CSTF3)

Product Code	CSB-BP618643HU
Abbreviation	CSTF3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12996
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SGDGATEQA AEYVPEKVKK AEKKLEENPY DLDAWSILIR EAQNQPIDKA RKTYERLVAQ FPSSGRFWKL YIEAEIKAKN YDKVEKLFQR CLMKVLHIDL WKCYLSYVRE TKGKLPYSYKE KMAQAYDFAL DKIGMEIMSY QIWVDYINFL KGVEAVGSYA ENQRITAVRR VYQRGCVNPM INIEQLWRDY NKYEEGINIH LAKKMIEDRS RDYMNARRVA KEYETVMKGL DRNAPSVPPQ NTPQEAQQVD MWKKYIQWEK SNPLRTEDQT LITKRVMFAY EQCLLVLGHH PDIWYEEAAQY LEQSSKLLAE KGDMNNAKLF SDEAANIYER AISTLLKKNM LLYFAYADYE ESRMKYEKVH SIYNRLLAIE DIDPTLVYIQ YMKFARRAEG IKSGRMIFKK AREDTRTRHH VYVTAALMEY YCSKDKSVAF KIFELGLKKY GDIPEYVLAY IDYLSHLNED NNTRVLFERV LTSGSLPPEK SGEIWARFLA FESNIGDLAS ILKVEKRRFT AFKEEYEGKE TALLVDRYKF MDLYPCSASE LKALGYKDVS RAKLAAIIPD PVPVAPSIVPV LKDEVDKPE YPKPDTQQMI PFQPRHLAPP GLHPVPGGVF PVPPAAVVLML KLLPPPICFQ GPFVQVDELM EIFRRCKIPN TVEEAVRIIT GGAPELAVEG NGPVESNAVL TKAVKRPNE SDEDEEKGA VPPVHDIYRA RQQKRIR
Source	Baculovirus
Target Names	CSTF3
Protein Names	Recommended name: Cleavage stimulation factor subunit 3 Alternative name(s): CF-1 77 kDa subunit Cleavage stimulation factor 77 kDa subunit Short name= CSTF 77 kDa subunit Short name= CstF-77
Expression Region	2-717
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.