



Recombinant Human Four and a half LIM domains protein 3 (FHL3)

Product Code	CSB-EP618796HU-B
Abbreviation	FHL3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13643
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	SESFDCAKC NESLYGRKYI QTDSGPYCVP CYDNTEFANTC AECQQLIGHD SRELFYEDRH FHEGCFRCCR CQRSLADEPF TCQDSELLCN DCYC SAFSSQ CSACGETVMP GSRKLEYGGQ TWHEHCFLCS GCEQPLGSR FVPDKGAHYC VPCYENKFAP RCARCSKTLT QGGVTYRDQP WHRECLVCTG CQTPLAGQQF TSRDEDPYCV ACFGELFAPK CSSCKRPIVG LGGGKYVSFE DRHWHHNCFS CARCSTSLVG QGFVPDGDQV LCQGCSQAGP
Source	E.coli
Target Names	FHL3
Protein Names	Recommended name: Four and a half LIM domains protein 3 Short name= FHL-3 Alternative name(s): Skeletal muscle LIM-protein 2 Short name= SLIM-2
Expression Region	2-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.