



Recombinant Human Activating signal cointegrator 1 (TRIP4)

Product Code	CSB-BP618991HU
Abbreviation	TRIP4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q15650
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	AVAGAVSGE PLVHWCTQQL RKTFLDVSE EIIQYVLSIE SAAEIREYVT DLLQGNQEGKK GQFIEELITK WQKNDQELIS DPLQQCFKKD EILDGQKSGD HLKRGRKKGR NRQEVPAFTE PDTTAEVKTP FDLAKAQENS NSVKKKTKFV NLYTREGQDR LAVLLPGRHP CDCLGQKHKL INNCLICGRI VCEQEGSGPC LFCGTLVCTH EEQDILQRDS NKSQKLLKLL MSGVENSQKLV DISTKDLLPH QELRIKSGLE KAIKHKDKLL EFDRTSIRRT QVIDDESDFY ASDSNQWLSK LERETLQKRE EELRELHAS RLSKKVTIDF AGRKILEEEN SLAEYHSRLD ETIQAIANGT LNQPLTKLDR SSEEPLGVLV NPNMYQSPQ WVDHTGAASQ KKAFRSSGFG LEFNSFQHQL RIQDQEFQEG FDGGWCLSVH QPWASLLVRG IKRVEGRSWY TPHRGRWIA ATAKKPSPE VSELQATYRL LRGKDVDFPN DYPGCLLGC VDLIDCLSQK QFKEQFPDIS QESDSPFVFI CKNPQEMVVK FPIKGNPKIW KLDSKIHQGA KKGLMKQNKV V
Source	Baculovirus
Target Names	TRIP4
Protein Names	Recommended name: Activating signal cointegrator 1 Short name= ASC-1 Alternative name(s): Thyroid receptor-interacting protein 4 Short name= TR- interacting protein 4 Short name= TRIP-4
Expression Region	2-581
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.