



# Recombinant Human Forkhead box protein D1 (FOXD1)

<b>Product Code</b>	CSB-YP619085HU
<b>Abbreviation</b>	FOXD1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q16676
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTLSTEMSDA SGLAEETDID VVGEGEDEED EEEEDDDEGG GGGPRLAVPA QRRRRRRSYA GEDELEDLEE EEDDDDILLA PPAGGSPAPP GPAPAAGAGA GGGGGGGGAG GGSAGSGAK NPLVKPPYSY IALITMAILQ SPKKRLTLSE ICEFISGRFP YYREKFPWQ NSIRHNLSLN DCFVKIPREP GNPGKGNYSWT LDPESADMFD NGSFLRRRKR FKRQPLLPPN AAAAESLLLR GAGAAGGAGD PAAAAALFPP APPPPPHAYG YGPYGCGYGL QLPPYAPPSA LFAAAAAAAA AAAFHPSPP PPPPPHGAAA ELARTAFGYR PHPLGAALPG PLPASAAKAG GPGASALARS PFSIESIIGG SLGPAAAAAA AAQAAAAAQA SPSPSPVAAP PAPGSSGGGC AAQAAVGPAALTRSLVAAA AAAASSVSSS AALGTLHQGT ALSSVENFTA RISNC
<b>Source</b>	Yeast
<b>Target Names</b>	FOXD1
<b>Protein Names</b>	Recommended name: Forkhead box protein D1 Alternative name(s): Forkhead-related protein FKHL8 Forkhead-related transcription factor 4 Short name= FREAC-4
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This intronless gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. The specific function of this gene has not yet been determined; however, it may play a role in tumor formation.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.