



Recombinant CCR4-NOT transcription complex subunit 7 (ccf-1)

Product Code	CSB-BP619134CXY
Abbreviation	ccf-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17345
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MASSSSGGAG GAGGASGAPE VKIHNVMYSN VEEEFARIRG FVEDYPYVAM DTEFPGVVAT PLGTFRSKED FNYQQVFCNV NMLKLIQVGF AMVNDKGELP PTGDVWQFNF NFSFAEDMFS HESVEMLRQA GIDFTLLQNN GIPTAVFGEL LTTSGELITDP RITWLTFSFG YDFGYLLKSI TLGDLPKEES TFFMCHKTLF PTSFDIKILL RTPNCASAKL KGGLQEVDQ LDVQRQGVHR QAGSDALLTA ATFFKIKKQF FGDNWNQIAP LICGHMFGLG SSSLFHHSSG STSRLGDETP QGLIGVPQQA
Source	Baculovirus
Target Names	ccf-1
Protein Names	Recommended name: CCR4-NOT transcription complex subunit 7 Alternative name(s): CCR4-associated factor 1 Short name= CAF1
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.