



# Recombinant Mitogen-activated protein kinase pmk-1 (pmk-1)

<b>Product Code</b>	CSB-MP619138CXY
<b>Abbreviation</b>	pmk-1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q17446
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MFPQTTMDHI LHPTPREGYY VVELNRSVWV VPNYIINLTP IGTGAYGTVC AAECTRSRTR VAIKKFNRPF QSIIHARRTY RELRLLRCMC HENIIDLLDV FTPNENVNDI EDVYFVSMML GADLSNILKI QRLNDDHIQF LUYQILRGLK YIHSADIIHR DLKPSNIAVN EDCELKILDF GLARQTDSEM TGYVATRWYR APEIMLNWMH YTQTVDVWSV GCILAEITG KTLFPGSDHI DQLTRIMSVT GTPDEEFLKK ISSEARNYI RNLPKMTRRD FKRLFAQATP QAIDLLEKML HLDPDRRPTA KEAMEHEYLA AYHDETDEPI AEEMDLNDDV RADTIDEWKK IIWEEISDFQ KNVAFADEEE DEEKMES
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pmk-1
<b>Protein Names</b>	Recommended name: Mitogen-activated protein kinase pmk-1 EC= 2.7.11.24 Alternative name(s): Stress-activated protein kinase pmk-1 p38 MAP kinase 1
<b>Expression Region</b>	1-377
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.