



Recombinant Deoxyribonuclease-2 (nuc-1)

Product Code	CSB-EP619142CXY
Abbreviation	nuc-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17778
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	AFSCKDQSG NDVDWFAVYK MPIEKDDGSV TGLAGGVAWY YVDVNKKGTL TPSAKTLDDN DQAIAYTLQQ YYDKQNDKTI FHVMYNDEPW GSKSTSGIKL EEILSNRVYS NYTHEDDSTS TAFGHTKGTI FFDGTSGVWL VHSVPLFPNP TKYEYPVSGH DYGQTMLCMT FKYAQLKSIG TQLFFNRPN IYSSNLPTNMA ADNADLAKAI AGQYQKGQPF QSVIELETMA GYSFTNFAKS KEFNADLYDT LVAPTLKTDL VVETWRRGSE IPLDCKLTYH ANDALSIHVG STTAFSYTKD HSKMAHSADM TKPWVCIGDI NRMTSQYVRG GGTTCISSSF LWKAYSVIAT QNNCA
Source	E.coli
Target Names	nuc-1
Protein Names	Recommended name: Deoxyribonuclease-2 EC= 3.1.22.1 Alternative name(s): Deoxyribonuclease II Short name= DNase II
Expression Region	22-375
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.