



Recombinant Human Exocyst complex component 3-like protein 4 (EXOC3L4)

Product Code	CSB-BP619168HU
Abbreviation	EXOC3L4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q17RC7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MPSPQTDTPG PELQSPKEAE EPQTPAQGSR RTSSRKEPNA HRKDGTRLGL GSLRQAFSRA SQRALTQVSK EDTGLFRRSS CSLFRSFRQA LNDGPATGHS QATPEVPSGV MNGVSQCAST GAASEELKPE AEGKSVADLI TERQLLAAFE QLLRLETLV AEKASRTFEQ DPTAFARRAM DVCLLYDGLA AEIGAIVRET LDSDGVDAVA LAELARVVA EEEAHPSPPD DGDFLRTPRR WRQHWEEAVR RSAQERVRRP GAGWAFGEAE GASGLAQLLA ELGGLVRRDL QKVRQEVQPA YAAAGFPAWE VYLRAFHSAV AQRLQELARD ARGCEQLYIL LDWAANVYGS PDLGAPGLA LPAEPLPPLL APDVARLES DYTSFLEAKI ASCFDSILQL EQSHWAAAEV PEVLQGLYQA PLSMDVHMLV AEHVKAAGAI SAELEATTLR ICTRALGLFV PRFEKAFLAS EAVSEPHLGA YINACEELRT SLLSRFPGTQ EELEKPLVTA TCSFQKHLQ GLQRELQPLF RVVCTRDWLT QDWLHPLMDK VVTFAGHLQR VARPRAQETL QEVHRFVVRE YLARALRPRE RFRGMERMHG SQKMSLDAQA ISDTFQGLGS EATWLDQAIQ CVAEILGETY KDDIQRHLET LIRSYPDIRR DHILAILALR RLGRQRNQLH LQHTQDLLRA AAGAAGAEAP RGRVLFEEIK VPSAMAVLIT CV
Source	Baculovirus
Target Names	EXOC3L4
Protein Names	Recommended name: Exocyst complex component 3-like protein 4
Expression Region	1-722
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.