



# Recombinant Human Exocyst complex component 3-like protein 4 (EXOC3L4)

<b>Product Code</b>	CSB-EP619168HU-B
<b>Abbreviation</b>	EXOC3L4
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q17RC7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MPSPQTDTPG PELQSPKEAE EPQTPAQGSR RTSSRKEPNA HRKDGTRLGL GSLRQAFSRA SQRALTQVSK EDTGLFRRSS CSLFRSFRQA LNDGPATGHS QATPEVPSGV MNGVVSQAST GAASEELKPE AEGKSVADLI TERQLLAAFE QLLRLETLV AEKASRTFEQ DPTAFARRAM DVCLLYDGLA AEIGAIVRET LDSDGVDAAA LAELARVSA EEEAHPSPPD DGDFLRTPRR WRQHWEEAVR RSAQERVRRP GAGWAFGEAE GASGLAQLLA ELGGLVRRDL QKVRQEVQPA YAAAGFPAWE VYLRAFHSAV AQRLQELARD ARGCEQLYIL LDWAANVYGS PDLGAPGLA LPAEPLPPLL APDVARLES DYTSFLEAKI ASCFDSILQL EQSHWAAA EVPEVLQGLYQA PLSMDVHMLV AEHVKAAGAI SAELEATTLR ICTRALGLFV PRFEKAFLAS EAVSEPHLGA YINACEELRT SLLSRFPGTQ EELEKPLVTA TCSFQKHLQ GLQRELQPLF RVVCTRDWLT QDWLHPLMDK VVTFAGHLQR VARPRAQETL QEVHRFVVRE YLARALRPRE RFRGMERMHG SQKMSLDAQA ISDTFQGLGS EATWLDQAIQ CVAEILGETY KDDIQRHLET LIRSYPDIRR DHILAILALR RLGQRNQHL LQHTQDLLRA AAGAAGAEAP RGRVLFEEIK VPSAMAVLIT CV
<b>Source</b>	E.coli
<b>Target Names</b>	EXOC3L4
<b>Protein Names</b>	Recommended name: Exocyst complex component 3-like protein 4
<b>Expression Region</b>	1-722
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.