



Recombinant Uncharacterized F-box protein C27F2.7 (C27F2.7)

Product Code	CSB-EP619213CXY-B
Abbreviation	C27F2.7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q18262
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MAVQQKRRGL SFIDAKASNV KKIIVMASTV AFEFNKATLR LRRSLHISPR SSPEVQRKAT AGENSEVGSP ESSLSTSRCS KRDRKLCAEL SSFALYPIFG ALSPSHFSSP NLLFQSEKNS SRRSFEFMQL PDTDICQIMS FLDAQSLLNL SQTCSHLRQL CLAHEDNAGK RDVTSHEITI SFNQIHR RTE VRLKRETR TDRQFISNNI RGVLPFSRA LTKITFETTV FVTDWLDEIL QLHRENRLIP LSLLFTGGAL TKGHQVTGAD LRSITETEFV DFVTKLQPHL QEVQLSTSR I FKISTDPSQL LGMISMLSSF GIVYERPLH FYHEEISNAI ALWKSDPLSR CCDVYMRPH DVSSSEWIKL AGSIDESRID PYTDEV LVSQ ITIKHSFLVH IDLVFHFH
Source	E.coli
Target Names	C27F2.7
Protein Names	Recommended name: Uncharacterized F-box protein C27F2.7
Expression Region	1-408
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.