



# Recombinant Mouse Lactadherin (Mfge8)

<b>Product Code</b>	CSB-BP013752MO
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P21956
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASGDFCDS SLCLNGGTCL TGQDNDIYCL CPEGFTGLVC NETERGPCSP NPCYNDAKCL VTLDTQRGDI FTEYICQCPV GYSGIHCETE TNYYNLDGEY MFTTAVPNTA VPTPAPTPDL SNNLASRCST QLGMEGGAIA DSQISASSVY MGFMGLQRWG PELARLYRTG IVNAWTASNY DSKPWIQVNL LRKMRVSGVM TQGASRAGRA EYLKTFKVAY SLDGRKFEFI QDESGGDKEF LGNLDNNSLK VNMFNPTLEA QYIKLYPVSC HRGCTLRFEL LGCELHGCSE PLGLKNNTIP DSQMSASSSY KTWNLRAFGW YPHLGRLDNQ GKINAWTAQS NSAKEWLQVD LGTQRQVTGI ITQGARDFGH IQYVASYKVA HSDDGVQWTV YEEQGSSKVF QGNLDNNSHK KNIFEKPFMA RYVRVLPVSW HNRITLRLEL LGC
<b>Source</b>	Baculovirus
<b>Target Names</b>	Mfge8
<b>Protein Names</b>	Recommended name: Lactadherin Alternative name(s): MFGM Milk fat globule-EGF factor 8 Short name= MFG-E8 SED1 Sperm surface protein SP47 Short name= MP47
<b>Expression Region</b>	23-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.