



Recombinant Mouse Lactadherin (Mfge8)

Product Code	CSB-EP013752MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21956
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	ASGDFCDS SLCLNGGTCL TGQDNDIYCL CPEGFTGLVC NETERGPCSP NPCYNDAKCL VTLDTQRGDI FTEYICQCPV GYSGIHCETE TNYYNLDGEY MFTTAVPNTA VPTPAPTPDL SNNLASRCST QLGMEGGAIA DSQISASSVY MGFMGLQRWG PELARLYRTG IVNAWTASNY DSKPWIQVNL LRKMRVSGVM TQGASRAGRA EYLKTFKVAY SLDGRKFEFI QDESGGDKEF LGNLDNNSLK VNMFNPTLEA QYIKLYPVSC HRGCTLRFEL LGCELHGCSE PLGLKNNTIP DSQMSASSY KTWNLRAFGW YPHLGRLDNQ GKINAWTAQS NSAKEWLQVD LGTQRQVTGI ITQGARDFGH IQYVASYKVA HSDDGVQWTV YEEQGSSKVF QGNLDNNSHK KNIFEKPFMA RYVRVLPVSW HNRITLRLEL LGC
Source	E.coli
Target Names	Mfge8
Protein Names	Recommended name: Lactadherin Alternative name(s): MFGM Milk fat globule-EGF factor 8 Short name= MFG-E8 SED1 Sperm surface protein SP47 Short name= MP47
Expression Region	23-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.