



# Recombinant 60 kDa chaperonin 2 (groL2)

<b>Product Code</b>	CSB-EP619320EJE
<b>Abbreviation</b>	groL2
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q19NJ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O1:K1 / APEC
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAAKDIRFGE DARARMVRGV NVLANAVKAT LGPKGRNVVL EKSFGAPTIT KDGVSVAKEI ELADKFENMG AQMVKEVASK TSDNAGDGTT TATVLAQALI REGMKAVAAG MNPMDLKRGI DKAVTSAVEE LKKISKPCST SKEIAQVGS SANSDTDIGE LIAKAMDKVG KEGVITVEEG SGTENELDVV EGMQFDRGYL SPYFINNPQS MQAELEDPFI LLHDKKISNV RDLLPILEGV AKAGKPLLIV AEDVEGEALA TLVVNTIRGI VKVCAVKAPG FGDRRKAMLE DMAILTGTV ISEEVGLSLE KATINDLGRA KKVQVSKENT TIIDGAGDTA DIEARIKQIK AQIEETSDY DREKLQERVA KLAGGVAVIK VGAATEVEMK EKKARVEDAL HATRAAVEEG IVPGGGVALI RAKAAIAELK GANEDQNHGI AIALRAMEAP LREIVTNAGD EPSVVLNRVA EGTGAFGYNA ANGEFGDMIE FGILDPTKVT RSALQNAASI AGLMITTEAM VAEAPKKEEP AAPGGGMGGM GGMDF
<b>Source</b>	E.coli
<b>Target Names</b>	groL2
<b>Protein Names</b>	Recommended name: 60 kDa chaperonin 2 Alternative name(s): GroEL protein 2 Protein Cpn60 2
<b>Expression Region</b>	1-545
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.