



# Recombinant Ceratotherium simum Forkhead box protein G1 (FOXG1)

<b>Product Code</b>	CSB-EP619330DRG
<b>Abbreviation</b>	FOXG1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q1A1A1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Ceratotherium simum (White rhinoceros) (Square-lipped rhinoceros)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLDMGDRKEV KMIPKSSFSI NSLVPEAVQS DNHHASHGHH NSHHPQH HHHHHHHPPP PAPQPPPPP PPQQPPAPQ PSQARGVPAA DDDKGPQQLL LPPPPPPPA AALDGAKADG LGGKGEPGGG PGELAPVGPD EKEKGAGAGG EEKKGAGEGG KDGEKGKEGE KKNKYEKPP FSYNALIMMA IRQSPEKRLT LNGIYEFIMK NFPYYRENKQ GWQNSIRHNL SLNKCFVKVP RHYDDPGKGN YWMLDPSSDD VFIGGTTGKL RRRSTTSRAK LAFKRGARLT STGLTFMDRA GSLYWPMSPF LSLHHPRASS TLSYNGTTSA YPSHPMPYSS VLTQNSLGN HSFSTANGLS VDRLVNGEIP YATHHLTAAA LAASVPCGLS VPCSGTYSLN PCSVNLLAGQ TSYFFPHVPH PSMTSQSSTS MSARAASSST SPQAPSTLPC ESLRPSLPSF TTGLSGGLSD YFTHQNQGSS SNPLIH
<b>Source</b>	E.coli
<b>Target Names</b>	FOXG1
<b>Protein Names</b>	Recommended name: Forkhead box protein G1 Short name= FoxG1
<b>Expression Region</b>	1-486
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.