



Recombinant *Saccharomyces cerevisiae* Golgi to ER traffic protein 4 (GET4)

Product Code	CSB-EP619563SVG-B
Abbreviation	GET4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12125
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVPAESNAVQ AKLAKTLQRF ENKIKAGDYY EAHQTLRTIA NRYVRSKSYE HAIELISQGA LSFLKAKQGG SGTDLIFYLL EVYDLAEVKV DDISVARLVR LIAELDPSEP NLKDVITGMN NWSIKFSEYK FGDPYLHNTI GSKLLEGDFV YEAERYFMLG THDSMIKYVD LLWDWLCQVD DIEDSTVAEF FSRLVFNLYL ISNISFAHES KDIFLERFIE KFHPKYEKID KNGYEIVFFE DYSDLNFLQL LLITCQTKDK SYFLNLKNHY LDFSQAYKSE LEFLGQEQYFN IVAPKQTNFL QDMMSGFLGG SK
Source	<i>E.coli</i>
Target Names	GET4
Protein Names	Recommended name: Golgi to ER traffic protein 4 Alternative name(s): Guided entry of tail-anchored proteins 4
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.