



Recombinant *Saccharomyces cerevisiae* Epsin-1 (ENT1)

Product Code	CSB-BP619610SVG
Abbreviation	ENT1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q12518
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MSKQFVRS AK NLVKGYSSTQ VLVRNATSND NHQVSKDSL I ELAEKSYDSA DFFEIMDMLD KRLNDKGKYW RHIAKALTVI DYLRFGSEN CVLWCRENLY IIKTLKEFRH EDDEGIDQQGQ IVRVKAKELT ALLSDDERLN EERNMNIKGR NRKGRRRRGT GRSDENDDDL QRAISASRLT AEEDERRRKQ DEDYETALQL SKEEEELKRL QDLQRMQQQQ GQQQLQQPMY YDIFGNPITP EEYAQFQLQQ QQQQQQQLQ QQPMYYDVFG NPITPEELAQ FQQQQQLQEQ QYLASMQQQQ QAMSNPFAK SEQSSSPKR NQLVAASSPQ QLQQKQKQEP LIQNRTGNQS MTDKYSKLINE LLATGTGIDT FGNVGEARIP AQHTKTGTFI NSQGTGYRQV SDDPNHNPFL NSQYTGLPST SVVPTQTGYG FGNQSQQQSQ NNGSNRNGYT LIDL
Source	Baculovirus
Target Names	ENT1
Protein Names	Recommended name: Epsin-1
Expression Region	1-454
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.