



# Recombinant *Saccharomyces cerevisiae* Epsin-1 (ENT1)

<b>Product Code</b>	CSB-EP619610SVG-B
<b>Abbreviation</b>	ENT1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12518
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSKQFVRS AK NLVKGYSSTQ VLVRNATSND NHQVSKDSL I ELAEKSYDSA DFFEIMDMLD KRLNDKGKYW RHI AKALTVI DY LIRFGSEN CVLWCRENLY IIKTLKEFRH EDDEGIDQQGQ IVRVKAKELT ALLSDDERLN EERNMNIKGR NRKGRRRRGT GRSDENDDDL QRAISASRLT AEEDERRRQ DEDYETALQL SKEEELKRL QDLQRMQQQQ GQQQLQQPMY YDIFGNPITP EEYAQFQLQQ QQQQQQQQLQ QQPMYYDVFG NPITPEELAQ FQQQQQLQEQ QYLASMQQQQ QAMSNPFAK SEQSSSPKR NQLVAASSPQ QLQQKQKQEP LIQNRTGNQS MTDKYSK LNE LLATGTGIDT FGNVGEARIP AQHTKTGTFI NSQGTGYRQV SDDPNHNPFL NSQYTGLPST SVVPTQTGYG FGNQSQQQSQ NNGSNRGYT LIDL
<b>Source</b>	E.coli
<b>Target Names</b>	ENT1
<b>Protein Names</b>	Recommended name: Epsin-1
<b>Expression Region</b>	1-454
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.