



# Recombinant Human Nuclear factor 1 A-type (NFIA)

<b>Product Code</b>	CSB-BP619635HU
<b>Abbreviation</b>	NFIA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q12857
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MYSPLCLTQD EFHPFIEALL PHVRAFAYTW FNLQARKRKY FKKHEKRMSK EEERAVKDEL LSEKPEVKQK WASRLAKLR KDIRPEYRED FVLTVTGKKP PCCVLSNPDQ KGKMRRIDCL RQADKVVRLD LVMVILFKGI PLESTDGERL VKSPQCSNPG LCVQPHHIGV SVKELDLYLA YFVHAADSSQ SESPSQPSDA DIKDQPENGH LGFQDSFVTS GVFSVTELVV VSQTPIAAGT GPNFSLSDLE SSSYYSMSPG AMRRSLPSTS STSSTKRLKS VEDEMDSPE EPFYTGQGRS PGSGSQSSGW HEVEPGMPSP TTLKSEKSG FSSPSPSQTSLGTAFTQHH RPVITGPRAS PHATPSTLHF PTSPHQQPG PYFSHPAIRY HPQETLKEFV QLVCPDAGQQ AGQVGFLNPN GSSQGVHNP FLPTPMLPPP PPPPMARPVP LPVPDTKPPT TSTEGGAASP TSPTYSTPST SPANRFVSVG PRDPSFVNIP QQTQSWYLG
<b>Source</b>	Baculovirus
<b>Target Names</b>	NFIA
<b>Protein Names</b>	Recommended name: Nuclear factor 1 A-type Short name= NF1-A Short name= Nuclear factor 1/A Alternative name(s): CCAAT-box-binding transcription factor Short name= CTF Nuclear factor I/A Short name= NF-I/A Short
<b>Expression Region</b>	1-509
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.