



# Recombinant Human Krueppel-like factor 1 (KLF1)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP619760HU                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Abbreviation</b>      | KLF1                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                              |
| <b>Uniprot No.</b>       | Q13351                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Sequence</b>          | MATAETALPS ISTLTALGPF PDTQDDFLKW WRSEEAQDMG PGPPDPTEPP<br>LHVKSEDQPG EEEDDERGAD ATWDLDLLLT NFSGPEPGGA PQTALAPSE<br>ASGAQYPPPP ETLGAYAGGP GLVAGLLGSE DHSGWVRPAL RARAPDAFVG<br>PALAPAPAPE PKALALQPVY PGPAGSSGG YFPRTGLSVP AASGAPYGLL<br>SGYPAMYAP QYQGHFQLFR GLQGPAPGPA TSPSFLSCLG PGTVGTGLGG<br>TAEDPGVIAE TAPSKRGRRS WARKRQAAHT CAHPGCGKSY TKSSHLKAHL<br>RTHTGKPYA CTWEGCGWRF ARSDELTRHY RKHTGQRPFR<br>CQLCPRAFSR SDHLALHMKR HL |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Target Names</b>      | KLF1                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Protein Names</b>     | Recommended name: Krueppel-like factor 1 Alternative name(s): Erythroid krueppel-like transcription factor Short name= EKLF                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 1-362                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                             |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Target Details</b>    | This gene encodes a hematopoietic-specific transcription factor that induces high-level expression of adult beta-globin and other erythroid genes. The zinc-finger protein binds to the DNA sequence CCACACCCT found in the beta hemoglobin promoter. Heterozygous loss-of-function mutations in this gene result in the dominant In(Lu) blood phenotype.                                                                       |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                     |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients,                                                                                                                                                                                                                                                                                                                                                   |



storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.