



Recombinant Human Cullin-4A (CUL4A)

Product Code	CSB-BP619782HU
Abbreviation	CUL4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13619
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MADEAPRKGS FSALVGR TNG LTKPAALAAA PAKPGGAGGS KKLVIKNFRD RPRLPDNYTQ DTWRKLHEAV RAVQSSTSIR YNLEELYQAV ENLCSHKVSP MLYKQLRQAC EDHVQAQILP FREDSLDSVL FLKKINTCWQ DHCRQMIMIR SIFLFLDRTY VLQNSTLPSI WDMGLELFR T HIISDKMVQS KTIDGILLI ERERSGEAVD RLLRSLLGM LSDLQVYKDS FELKFLEETN CLYAAEGQRL MQEREVPEYL NHVSKRLEEE GDRVITYLDH STQKPLIACV EKQLLGEHLT AILQKGLDHL LDENRVPDLA QMYQLFSRVR GGQQALLQHW SEYIKTFGTA IVINPEKDKD MVQDLLDFKD KVDHVIEVCF QKNERFVNLM KESFETFINK RPNKPAELIA KHVDSKLRAG NKEATDEELE RTLDKIMILF RFIHGKDVFE AFYKKDLAKR LLVGKSASVD AEKSM LSKLK HECGAAFTSK LEGMFKDMEL SKDIMVHFKQ HMQNSQSDSGP IDLTVNILTM GYWPTYTPME VHLTPEMIKL QEVFKAFYLG KHSGRKLQWQ TTLGHAVLKA EFKEGKKEFQ VSLFQTLVLL MFNEGDGFSF EEIKMATGIE DSELRR TLQS LACGKARVLI KSPKGKEVED GDKFIFNGEF KHKLFRIKIN QIQMKETVEE QVSTTERVFQ DRQYQIDAAI VRIMKMRKTL GHNLLVSELY NQLKFPVKPG DLKKRIESLI DRDYMERDKD NPNQYHYVA
Source	Baculovirus
Target Names	CUL4A
Protein Names	Recommended name: Cullin-4A Short name= CUL-4A
Expression Region	1-759
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.