



Recombinant Human Four and a half LIM domains protein 1 (FHL1)

Product Code	CSB-YP619784HU
Abbreviation	FHL1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13642
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	AEKFDCHYC RDPLQGKKYV QKDGHHCCCLK CFDKFCANTC VECRKPIGAD SKEVHYKNRF WHDTCFRCAK CLHPLANETF VAKDNKILCN KCTTREDSPK CKGCFKAIVA GDQNVVEYKGT VWHKDCFTCS NCKQVIGTGS FFPKGEDFYC VTCHETKFAK HCVKCNKAIT SGGITYQDQP WHADCFVCVT CSKLAGQRF TAVEDQYYCV DCYKNFVAKK CAGCKNPITG KRTVSRVSHV VSKARKPPVC HGKRLPLTLF PSANLRGRHP GGERTCPSWV VVLYRKNRSL AAPRGPLVK APVWWPMKDN PGTTTASTAK NAP
Source	Yeast
Target Names	FHL1
Protein Names	Recommended name: Four and a half LIM domains protein 1 Short name= FHL-1 Alternative name(s): Skeletal muscle LIM-protein 1 Short name= SLIM Short name= SLIM-1
Expression Region	2-323
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.