



Recombinant Human RNA-binding protein 39 (RBM39)

Product Code	CSB-BP619875HU
Abbreviation	RBM39
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q14498
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	ADDIDIEAM LEAPYKKDEN KLSSANGHEE RSKKRKKS KS RSRSHERKRS KSKERKRSRD RERKKS KSRE RKR SR SKERR RSR SRSDRR FRGRYRSPYS GPKFNSAIRG KIGLPHSIKL SRRRSR SKSP FRKDKSPVRE PIDNLTPEER DARTVFCMQL AARIRPRDLE EFFSTVGKVR DVRMISDRNS RRSKGIAYVE FVDVSSVPLA IGLTGQRVLG VPIIVQASQA EKNRAAAMAN NLQKGSAGPM RLYVGS LHFN ITEDMLRGIF EPFGRIESIQ LMMDSETGRS KGYGFITFSD SECAKKALEQ LNGFELAGRP MKVGHV TERT DASSASSFLD SDELERTGID LGTTGRLQLM ARLAEGTGLQ IPPAAQQALQ MSGSLAFGAV AEFSFVIDLQ TRLSQQTEAS ALAAAASVQP LATQCFQLSN MFNPQTEEEV GWDTEIKDDV IEECNKHGGV IHIVVDKNSA QGNVYVKCPS IAAAIAAVNA LHGRWFAGKM ITAAYVPLPT YHNLFPDSMT ATQLLVPSRR
Source	Baculovirus
Target Names	RBM39
Protein Names	Recommended name: RNA-binding protein 39 Alternative name(s): Hepatocellular carcinoma protein 1 RNA-binding motif protein 39 RNA-binding region-containing protein 2 Splicing factor HCC1
Expression Region	2-530
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.