



# Recombinant Human RNA-binding protein 39 (RBM39)

<b>Product Code</b>	CSB-EP619875HU
<b>Abbreviation</b>	RBM39
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q14498
<b>Storage Buffer</b>	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ADDIDIEAM LEAPYKKDEN KLSSANGHEE RSKKRKKS KS RSR SHERKRS KSKERKRSRD RERKKS KSRE RKRSRSKERR RSRRSRDRR FRGRYRSPYS GPKFNSAIRG KIGLPHSIKL SRRRSRSKSP FRKDKSPVRE PIDNLTPEER DARTVFCMQL AARIRPRDLE EFFSTVGKVR DVRMISDRNS RRSKGIAYVE FVDVSSVPLA IGLTGQRVLG VPIIVQASQA EKNRAAAMAN NLQKGSAGPM RLYVGS LHFN ITEDMLRGIF EPFGRIESIQ LMMDSETGRS KGYGFITFSD SECAKKALEQ LNGFELAGRP MKVGHV TERT DASSASSFLD SDELERTGID LGTTGRLQLM ARLAEGTGLQ IPPAAQQALQ MSGSLAFGAV AEFSFVIDLQ TRLSQQTEAS ALAAAASVQP LATQCFQLSN MFNPQTEEEV GWDTEIKDDV IEECNKHGGV IHIVYDKNSA QGNVYVKCPS IAAAIAAVNA LHGRWFAGKM ITAAYVPLPT YHNLFPDSMT ATQLLVPSRR
<b>Source</b>	E.coli
<b>Target Names</b>	RBM39
<b>Protein Names</b>	Recommended name: RNA-binding protein 39 Alternative name(s): Hepatocellular carcinoma protein 1 RNA-binding motif protein 39 RNA-binding region-containing protein 2 Splicing factor HCC1
<b>Expression Region</b>	2-530
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

---

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.