



# Recombinant Human TNFAIP3-interacting protein 1 (TNIP1)

<b>Product Code</b>	CSB-MP619959HU
<b>Abbreviation</b>	TNIP1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q15025
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MEGRGPYRIY DPGGSVPSGE ASAAFERLVK ENSRLKEKMQ GIKMLGELLE ESQMEATRLR QKAEELVKDN ELLPPPSPSL GSFDP LAELT GKDSNV TASP TAPACPSDKP APVQKPPSSG TSSEFEVVTP EEQNSPESSS HANAMALGPL PREDGNLMLH LQRLETTLSV CAEEDPHGQL FTHLGRMALE FNRLASKVHK NEQRTSILQT LCEQLRKENE ALKAKLDKGL EQRDQAAERL REENLELKKL LMSNGNKEGA SGRPGSPKME GTGKKAVAGQ QQASVTAGKV PEVVALGAAE KVKMLEQQR SELLEVNKQW DQHFRRSMKQQ YEQKITELRQ KLADLQKQVT DLEAEREQKQ RDFDRKLLLA KSKIEMEETD KEQLTAEAKE LRQVKYLQD QLSPLTRQRE YQEKEIQR LN KALEEALSIQ TPPSSPPTAF GSPEGAGALL RKQELVTQNE LLKQQVKIFE EDFQRERSDR ERMNEEKEEL KKQVEKLQAQ VTLSNAQLKA FKDEEKAREA LRQQKRKAKA SGERYHVEPH PEHLGAYPY AYPPMPAMVP HHGFEDWSQI RYPPPPMAME HPPPLPNSRL FHLPEYTWRL PCGGVRNPNQ SSQVMDPPTA RPTEPESPKN DREGPQ
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TNIP1
<b>Protein Names</b>	Recommended name: TNFAIP3-interacting protein 1 Alternative name(s): HIV-1 Nef-interacting protein Nef-associated factor 1 Short name= Naf1 Nip40-1 Virion-associated nuclear shuttling protein Short name= VAN Short nam
<b>Expression Region</b>	1-636
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.